

# Cascade CMXErg+

Independently adjustable crank arms  
Bi-directional resistance  
Virtually zero watt starting resistance  
Handles can be interchanged with pedals for lower extremity training  
34 pound total weight- with built in carry handle  
Full commercial warranty- light commercial pricing!

CMX SMARTTECH POWER CONSOLE displays key cycling feedback:

- Level/Resistance 1-8
- Calories
- Time Exercised
- Speed - MPH or KPH
- Distance



Pedals come standard for lower extremity exercise.

Length: 16" (40.6 cm)  
Height: 18" (48.3 cm)  
Width: 13" (33 cm)  
Weight: 34 lbs. (15.4 kg)  
Shipping Weight: 38 lbs. (17.2 kg)



Medical Fitness Solutions

543 E. Alvarado St,

Fallbrook, CA 92028 800.831.7665