

THE ROAD TO RECOVERY STARTS HERE

L7 & L8 REHAB TREADMILLS

Made with the finest materials, the L7 and L8 Rehabilitation Treadmills provide the best and safest support for all your rehabilitative programs. With a low step-up height, zero starting speed, 300 micro-amp isolation leakage kit and extended parallel medical handrails, these treadmills are ideal for any physical therapy and sports conditioning center. A number of unique options are also available, such as a reversing treadbelt for downhill walking, 100 micro-amp hospital isolation leakage kit, and an orthopedic suspension system. For a complete list of options please visit us at www.landice.com.

TREADMILL FRAME COMPARISON

Solid aluminum construction – rust free and durable

	L8	L7
User Capacity	500 lbs.	400 lbs.
Equipment Weight	300 lbs.	340 lbs.
Dimensions	35" W x 83" L x 61.5" H	35" W x 76" L x 60" H
Running Surface	22" W x 63" L	20" W x 58" L
Power	110 VAC, 60Hz, 15A 220 VAC, 60Hz, 15A (optional)	110 VAC, 60Hz, 15A 220 VAC, 60Hz, 15A (optional)
Rollers	3 ½" 22 lbs.	2 ½" 14 lbs.
Step-up Height	7 ½"	5 ½"
Speed	0.1 to 12 mph	0.1 to 12 mph
Incline	0% to 15%	0% to 15%
Drive System	4 HP Continuous Duty Motor	4 HP Continuous Duty Motor
Deck	1" - Reversible	1" - Reversible
Control Panel	Rehabilitation	Rehabilitation
Heart Rate	Interactive Contact Heart Rate	Interactive Contact Heart Rate
Warranty*	5-year parts, 3-year console, 2-year wear items, 1-year labor	5-year parts, 3-year console, 2-year wear items, 1-year labor

*Treadmills are covered by the LTD warranty which offers coverage for usage of less than 5 hours per day.



Landice engineers, assembles and tests the finest specialty fitness equipment on the market and is one of the only manufacturers to do so in the USA. Landice has been recognized repeatedly by leading consumer rating publications as having the world's finest treadmills and ellipticals.