

# koreHbot™

Hyperbaric Chamber Lying Type

Model: ST801 / ST901



## Features

- 1.3 ATA / 1.4 ATA / 1.5 ATA operating pressure.
- 2 zipper seal for easier entrance.
- 7 large transplant viewing windows to prevent claustrophobia
- Easy to wash cotton chamber protection cover to maintain cleanliness.
- Most popular model for home treatment or commercial usage.
- Emergency valve - to speed up the depressurization in emergency.
- Easy operation - one person can operate without assistance.
- Delivers 93% oxygen unhdre pressure via nasal cannula or oxygen mask.
- Internal pressure gauge - user can observe the pressure while inside the chamber
- Internal metal frame - holds the shape when deflated

Model		ST801	ST901
Chamber Capsule	Pressure	1.5 ATA	1.4 ATA
	Material	TPU	TPU
	Size (DxL)	80x225cm (32x89in)	90x225cm (36x89in)
	Weight	13kg (lbs)	16kg (lbs)
Air Compressor	Size	39x24x26cm (16x10x11in)	
	Type	Quiet & Oil Free type. Two super absorption activated carbon filters in air input and outlet to protect.	
Oxygen Concentrator	Size	5 Liter/min / 10 Liter/min	
	Flow Rate	5 Liter/min / 10 Liter/min	
	Description	PSA molecular sieve high technology. Continuous oxygen production, no need for oxygen tank.	
Air Cooler	Size	18x12x35cm (7.2x5x14in)	
	Description	Cooling and dehumidification	