



- Powered by Electricity - No Liquid Nitrogen!
- Installation & Delivery Included - Seamless & Professional
- Extreme Low Temperatures
- Fully Automated
- Communication Capability
- Pre-set On/Off
- Built-In Analytics
- Designed, Built and Supported in San Diego - Visit our Showroom!

Features

- **Temperature Range:** -100°F to -115°F, with the capability to reach as low as -140°F/-180°F
- **Cooling Time:** Approx. 45-60 minutes
- **Safety Features:** Temperature Sensors & Safety Interlocks
- **Utility Requirements:** 208-230 VAC, 3-ph, 40 Amp, 5-wire
- **Product Support:** Smart device access capability
- **Dimensions Required: Cryochamber**
 - **Module Space:** 54 inches (1372 mm) x 63 inches (1600.2 mm) x 94 inches (2387.6 mm)
 - **Operating Space:** 78 inches (1982 mm) x 87 inches (2210.2 mm) x 118 inches (2997.6 mm)
 - **Approximate Installation Space:** 102 inches (2592 mm) x 111 inches (2820.2 mm) x 142 inches (3607.6 mm)