


Korebalance Active

Model Number: Korebalance Active



*Laptop Stand, Support Rails, and Digital Pump Shown are Optional



We're driven because you're driven. The Korebalance™ Active is designed for use by medical centers and clinics, medical and PT offices, assisted living facilities, recreation centers, health and fitness clubs, universities and professional sport teams.

A Study in Perfection. The Korebalance™ System uses the latest in virtual and interactive technology, offering innovative high-tech balance assessment and training. Backed by over two decades of medical experience, the Korebalance™ systems will provide you and your patients with the best possible documented neurosensory assessment and training.

SPECIAL FEATURES:

- 17" Touchscreen Computer
- Optional:
 - Adjustable Laptop Stand
 - Carry Case
 - Support Rails
 - Rechargeable Digital Pump
- Variable Stability Platform with 23" (58.4 cm) Diameter
- Includes Korebalance® testing/training software suite & games.
- Optional Cognitive, Parkinson's Software Suite
- Patient Weight Capacity: 500 lbs (227 kg)
- Product Weight: 48 lbs (22 kg)
- Software Available in English, Spanish, French, German and Mandarin
- Shipping Weight Total: 56 lbs (25 kg)
- Shipping Dimensions: 27" x 27" x 9": 28" x 27" x 9"

Warranty: Frame: 5 year; Computer: 1 years; Parts: 1 year

In the interest of product development and improvement, specifications subject to change without notice.