



The 7.0T treadmill is engineered to deliver superior rehabilitation and patient results. With a comprehensive suite of features, users can enhance their therapy sessions like never before.

This treadmill provides exceptional flexibility, allowing you to tailor every session to meet each client's unique needs:

- Full-length, adjustable height handrails for added stability
- -10% - 15% decline/incline range
- Low starting speed and step-up block for easy on/off access
- 3.0 HP motor with bidirectional walking direction capabilities

The 7.0T is designed to withstand intense therapeutic sessions, making it the ideal choice for any commercial therapy setting.



Intuitive, easy-to-read display



Fully adjustable parallel bars support a wide variety of users



Front and rear motors allow for true up/down walking

## Specifications

|  |   |
|--|---|
| Incline  | -10% to 15%   |
| Bidirectional Speed  | Forward 0.1 - 10 mph and reverse 0.1 - 3 mph in 0.1 increments  |
| Motor  | 3HP continuous duty, IEC Class S6   |
| Belt and Deck  | 22" x 60" self-lubricating deck   |
| Power  | 15 volts AC, 60Hz, 20 amps, NEMA-20P  |
| Max User Weight  | 440 lb / 200 kg   |
| Dimensions   | 94.8" (83.7" w/o rear step) x 48.7" x 57.5"   |
| Product Weight   | 496 lb  |
| Compliance   | CAN/CAS-C22.2 No. 60601-1:14, ANS/AAMI ES60601-1:2005+A2 (R2012) +A1<br>IEC 60601-1-2:2014, EN 60601-1-2:2015<br>IEC 60601-1:2015+A1:2012, EN 60601-1:2006+A1:2013+A12:2014<br>EN ISO 20957:2013, MDD 93/42/EEC Class IIa, NB: 0123 |
| Warranty - Commercial (Dues paying)/<br>Rehab / Light Commercial | 5 years parts, 1 year labor   |