



MEDICAL MULTI HIP WITH ROM  
MFS-143-M-MH



**SPECIAL FEATURES:**

- 40 oz. seamless, antibacterial upholstery for superior durability, comfort, and easy cleaning
- Powder-coated frame with a scratch-resistant finish — available in black, white, silver, or custom colors
- Heavy-duty 11-gauge steel frame engineered for lasting strength and stability
- Heavily padded rollers provide superior comfort during exercise
- Assist handles with custom rubber grips for a secure, ergonomic hold
- Large non-slip foot platform ensures stability and safety
- Adjustable carriage accommodates a variety of user heights and training positions
- Adjustable leg pad arm allows precise positioning for optimal comfort and performance

In the interest of product development and improvement, specifications subject to change without notice.