


Korebalance Standard

Model Number: KB-Standard



We're driven because you're driven. The Korebalance™ Standard is designed for use by medical centers and clinics, medical and PT offices, assisted living facilities, recreation centers, health and fitness clubs, universities and professional sport teams.

A Study in Perfection. The Korebalance™ System uses the latest in virtual and interactive technology, offering innovative high-tech balance assessment and training. Backed by over two decades of medical experience, the Korebalance™ systems will provide you and your patients with the best possible documented neurosensory assessment and training.

SPECIAL FEATURES:

- 24" Touchscreen Computer
- Adjustable Monitor Mount Arm
- Fixed Safety Rails
- Platform has 360° horizontal & 20° vertical range of movement.
- Includes Korebalance® testing/training software suite & games.
- Optional Cognitive, Parkinson's Software Suite
- Exercise Band Attachment Points
- Patient Weight Capacity: 500 lbs (227 kg)
- Machine Weight: 315 lbs (143 kg)
- Footprint: 44"x58" (112 cm x 148 cm)
- Max Height: 80" (203 cm)
- The arms at the widest (at entry) is 35" and at the narrowest is 24"
- Power: 100-240 VAC; 50/60 Hz; 500 Watt
- Software Available in English, Spanish, French, German and Mandarin
- Frame Color: White
- Shipping Weight: 420 lbs (191 kg)
- Shipping Dimensions: 63" x 45" x 81"

Warranty: Frame: 15 year; Computer: 1 years; Parts: 1 year

In the interest of product development and improvement, specifications subject to change without notice.